



	Weight (g / bag)	Applicable Raw Material Weight (kg)	Material of the Bag
Koji Starter for Ginjo Sake	70	200	Nylon Plastic Bag
Koji Starter for Sake and Mirin	70	200	Nylon Plastic Bag
Koji Starter for Shochu	70	200	Nylon Plastic Bag
Koji Starter for Miso and Amazake	30	200	Nylon Plastic Bag
Koji Starter for Soy Sauce	40	240	Nylon Plastic Bag

These Products are Powder Type.

Only Koji spores are collected and mixed with extender (α -starch).

(Extender is added with the Powder Type Koji Starter to improve the spore engraftment rate)